

Real science with handheld LIBS


LIBZ Technology can see so many more elements compared to other handheld elemental technologies. From **hydrogen to uranium**. With Alloy, Geo, & Empirical Apps, you can confidently access and fine tune factory standard calibrations or develop your own custom profiles. With SciAps' eye-safe laser there are **no x-rays**, and no regulations to worry about!



Unleash the power of LIBZ handheld.

Create your own profiles with SciAps Profile Builder App. Run it on a PC, laptop or tablet. Locate elemental lines after testing samples, calculate intensities and intensity ratios, performs spectral overlays, baseline subtractions and smoothing. Generate calibration methods using reference samples with known assays.

Test unknown samples - Get a qualitative elemental breakdown, without calibrations. Want quantitation? Enter known assays and build calibration curves



SciAps.com
sales@sciaps.com
Woburn, MA USA
+1 339.927.9455

SciAps

LIBZ
TECHNOLOGY



PITTCONTM
CONFERENCE & EXPO 2015